


# INOGEN ONE<sup>®</sup> G5<sup>™</sup> System



Inogen One  
AC Power Supply\*  
Item #BA-501



Inogen One G5  
DC Power Cable\*  
Item #BA-306



Inogen One G5  
Lithium Ion Battery\*  
Up to 6.5 hours run time  
Item #BA-500



Inogen One G5  
Carry Bag\*  
Item #CA-500



Inogen One G5  
Backpack\*\*  
Item #CA-550



Inogen One<sup>®</sup> G5  
Extended Life Lithium Ion Battery\*\*  
Up to 13 hours run time  
Item #BA-516

## The Inogen One G5 System (IS-500):

- Compact and Lightweight
- Flow Settings 1-6
- Intelligent Pulse Oxygen Delivery
- Designed for 24 Hour Use
- Easy-to-read LCD display
- Inogen Connect App

\*Included with system

\*\*Accessory sold separately



## Oxygen. Anytime. Anywhere<sup>®</sup>

Once again, Inogen innovation sets a new standard in oxygen delivery.

The Inogen One G5 offers the most oxygen per pound for a portable oxygen concentrator on the market today. With flow settings from 1-6, the Inogen One G5 is designed to dramatically increase independence for most users 24/7. The Inogen One G5's user interface has been designed with an easy-to-read LCD display and simple control functions. It can be charged at home or on the go - giving you the freedom of Oxygen. Anytime. Anywhere.

The Inogen One G5 is your one solution, your single solution, for oxygen at home or away.

Weight	2.2 kg. (includes single battery)
Size	Length: 7.19" Width: 3.26" Height: 8.15" (includes single battery)
Oxygen Flow	Pulse dose delivery system Six flow settings 1, 2, 3, 4, 5, 6
Power	AC Power Supply 100-240V, 50-60Hz (auto sensing to allow worldwide use) DC Power Cable: for mobile use in car
Battery	Duration (single battery): Up to 6.5 hours Recharge up to 3 hours with AC or DC power Duration (double battery): Up to 13 hours Recharge up to 6 hours with AC or DC power
Noise Level	38 dBA*
Warranty	2 year warranty**
Operation	Simple control functions and easy-to-read LCD display Inogen Connect App
Use	Designed for 24/7 use and conforms to all applicable FAA requirements for POC use on board an aircraft

\*At flow setting 2

\*\*See Inogen One<sup>®</sup> Limited Warranty Statement for complete warranty information

Manufactured by:

  
326 Bollay Drive  
Goleta, California 93117  
[www.inogen.com](http://www.inogen.com)



## Marketed by

Office No. 506, 5th Floor, Office Tower – B, Elante,  
Plot no. 178-178A, Ind Area, Phase I, Chandigarh

Call : +91 7341103371  
[www.sanraishop.com](http://www.sanraishop.com)