

OXLIFE INDEPENDENCE[®]

YOUR OXYGEN PRESCRIPTION AT HOME OR OUTDOORS

- 12 pulse settings and 3 LPM continuous flow, giving you total flexibility to comfortably travel ***Best in Class***
- Swappable batteries allowing you to extend the battery time without needing to charge the entire device
- 2 in 1: Functions as a both a portable and stationary device in only 1 device
- Robust, reliable, and durable with a sturdy frame and inbuilt wheels
- Transportable: Easy to carry in a car, train, or plane
- Longest battery run time ***Best in Class***



SUPERIOR PRODUCT FEATURES

- Industry Leading Durability and Reliability = Lower Operating Costs
- Energy Smart Technology® = Longest Battery Life in Class
- Easy Battery Swap While Device is in Operation
- Lightweight Metal / Polycarbonate with ***Best in Class*** Components
- Integrated Cart with 6-inch Wheels Offers ***Best in Class*** Maneuverability

PRODUCT SPECIFICATIONS¹

Dimensions	Height: 20.2" • Width: 10.9" • Depth: 9.45"
Weight	7.5 kgs (with cart & battery)
Battery Recharge Time	2.5 hours
Oxygen Concentration	87% - 95% at sea level at all settings
DC Power	12 - 15 VDC • 15 amps
AC Power	Universal Input: 100 - 240 VAC 50 - 60 Hz • 2.0A maximum
Battery Duration	5.75 hours @ 2 pulse at 20 bp 2.5 hours @ 2 LPM continuous
Oxygen Flow Rate (Pulse 15-40 BPM)	Continuous: 0.5 - 3 LPM in 0.5 increments Pulse: 0.5 - 6.0 in 0.5 increments (8 - 96 ML)

Customisable: Additional batteries, car charger, and external battery charger are available for extended travel

